WINBACK

BETTER FASTER GENTLER

FR

EN

ES

IT

BACK

USER MANUAL

BACK 

WINBACK

BETTER FASTER GENTLER

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ENGLISH (EN)
USER MANUAL.....p.29

ESPAÑOL (ES)
MANUAL DE USARIO.....p.52

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BACK

REHABILITATION

SPORT

WOMEN

AESTHETIC

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WOMEN

AESTHETIC

Depending on the software updates for the BACK4 device, this manual might be subject to change.

To consult the latest version, scan the QR code opposite.

TABLE OF CONTENT

1. GENERAL INFORMATION	p.30
1.1 Introduction	p.30
1.2 Training.....	p.30
1.3 Warranty	p.30
2. SAFETY	p.31
2.1 Contraindications and side effect.....	p.31
2.1.1 Contraindications	p.31
2.1.2 Side effects.....	p.31
2.2 Precautions for use	p.31
3. DEVICE DESCRIPTION.....	p.32
3.1 Location of BACK4 controls	p.32
3.2 Main screen/functions	p.33
3.3 BACK4 energy	p.34
3.4 Accessories.....	p.35
4. OPERATING PRINCIPLES	p.37
4.1 Effects of the high-frequency current	p.37
4.2 Modes and functions.....	p.37
5. PERFORMING A TREATMENT	p.39
5.1 Preparing the device and performing a treatment.....	p.39
5.2 Treatment	p.39
6. MAINTENANCE OF THE DEVICE	p.42
6.1 General maintenance	p.42
6.2 Remote maintenance.....	p.43
7. TROUBLESHOOTING & ERROR MESSAGES	p.44
7.1 Troubleshooting	p.44
7.2 Error messages	p.45
8. TECHNICAL DATA.....	p.45
8.1 Technical specification	p.45
8.2 Composition conformity statement	p.48
8.3 Technical maintenance.....	p.48
9. END OF LIFE.....	p.48
10. SAFETY AND ELECTRICAL STANDARDS	p.48

1. GENERAL INFORMATION

This user manual contains essential safety information for the use of the BACK4 tecartherapy device and TECAR therapy. Please read and follow all safety instructions before using this device.

- Learn about the latest clinical developments and appropriate indications and contraindications.
- Observe all applicable precautionary measures for the treatment.

1.1 INTRODUCTION

BACK4 is the latest medical device from SWIMS S.A.S./WINBACK. This new device uses different characteristics of electrical currents with pulsed and modulated signals to create multiple therapeutic effects.

- TECAR (high frequency), stimulating the body's natural healing mechanisms and promoting cellular exchanges.
- Hi-TENS (high frequency with low-frequency pulses), increasing the analgesic power of TENS.
- Hi-EMS (medium frequency modulated in low frequencies), deep muscle contraction to obtain muscle strengthening and drainage.

BACK4 uses the patented TECAR^{core} and Hi-EMS^{core} technologies. It can be used to help patients recover and resume daily activities faster and in the best condition. It is suitable for use on patients aged 18 to 80, in good general health, presenting with an indication for which the device is suitable and without any contraindications linked to treatment.

The indications for the BACK4 medical device are as follows: relief of muscular and symptomatic pain (acute and chronic) - relaxation of muscle spasms - improvement of bone healing linked to osteoarthritis pathologies - improved joint function and mobility - management of muscular and tendon pathologies due to trauma - improvement of microcirculation and the healing process - muscle strengthening with increased flexibility.

BACK4 is designed for doctors and physiotherapists and is suitable for use in medical offices, clinics, hospitals and rehabilitation centers. It must be used by trained professionals with appropriate monitoring, and is for external use only. It must not be used in a home environment.

BACK4 is approved for use in all EU countries, sans exception.

1.2 TRAINING

Following delivery of your medical device, you will receive a one-day training session from a health professional appointed by the company SWIMS S.A.S./WINBACK. Do not hesitate to contact our service by phone on +33 (0)4 83 66 16 66 to arrange a date for your training.

1.3 WARRANTY

Your device is guaranteed from the date of delivery against any material defect inherent in its manufacture for a period of two years for a new device and one year for a second-hand device. Accessories delivered with the device are guaranteed for one year (with the exception of second-hand accessories, which are guaranteed for 6 months). SWIMS S.A.S./WINBACK may choose, at its discretion, to repair or replace defective parts during the warranty period, on the sole condition that the repairs and/or replacements are carried out by SWIMS S.A.S./WINBACK or one of its authorized representatives.

Note: The duration of the warranty may vary according to the regulation of each country.

SWIMS S.A.S./WINBACK declines all responsibility in the following cases:

- For the safety, reliability, and performance of the product if the set-up, settings, modifications, and/or any repairs are carried out by non-authorized persons.
- If the security label has been removed and/or moved.
- If the device is not used in accordance with the user manual.
- If the equipment is used without conductive cream and/or with a cream that is different from the one supplied by SWIMS S.A.S./WINBACK.
- If the installation or electrical network does not meet national or European requirements.
- Regarding normal wear and tear of accessories (plate, cables, electrodes).
- In the event of intentional or negligent damage to the equipment or its accessories.
- Penetration of liquids and creams into the equipment and/or the accessories.
- If the user does not follow the cleaning and maintenance protocols set out in this manual.
- Intentional misuse of the device (for example to treat a patient with contraindications; if the emergency button is unplugged during treatment; if the device is opened, etc.)

2. SAFETY

2.1 CONTRAINDICATIONS AND SIDE EFFECTS

2.1.1 Contraindications

The practitioner must refrain from any treatment for the following cases:

- Pregnant women
- Minor
- Cancer and cancerous lesions
- Skin disorders (Eczema, Shingles) - Skin burns - Open wounds
- Coagulation disorders – Phlebitis – Thrombophlebitis
- Insensitivity to heat
- Local insensitivity to pain and/or touch
- Fever – Bacterial infection – Infectious disease
- Severe hypertension or hypotension
- Electrical implants (pacemaker, insulin pump, neurostimulator)

The practitioner must refrain from any treatment involving direct application (direct contact with the electrode) on the following areas: the brain area, the eyes and around the eyes, and the heart area. We recommend the practitioner to make patients sign an informed consent before the first treatment. You should inform the patient about how the device works, the expected benefits, possible side effects, warnings, precautions and contraindications.

2.1.2 Side effects

- In some cases, a transient increase in pain may be felt within 24 hours after the session if the intensity is too high. This feeling fades spontaneously. In this case, it is advisable to work in athermal mode (low intensity) during the first three sessions.
- In very rare cases, the patient may be subject to an allergy to the conductive cream or superficial burns.
- In rare cases when the entire body is being treated, reactive hypotension may be observed. If the patient is lying down, be sure to raise them gradually after the treatment. In case of marked hypotension, hyperthermia is not recommended.

Important notes:

Prefer athermal treatments in the following cases:

- In case of treatment on a bone area (the thermal effect being faster).
- In case of anatomical vascular damage.
- In case of major inflammation, during the first three sessions.

In the event of any serious incident occurring in connection with the device, please inform the SWIMS S.A.S./WINBACK company by contacting the customer service (sav@swims.team). SWIMS S.A.S./WINBACK will then inform the manufacturer and the competent authority of the Member State in which the user and/or the patient is based.

2.2 PRECAUTIONS FOR USE

Operation:

- Do not operate the device when treating patients who are connected to other equipment.
- The accessories must always be entirely in contact with the patient's skin.
- Never use a sharp or pointed object like the tip of a pen or pencil on the screen.
- The accessories must be used with care. Inappropriate handling may damage them and impact their proper functioning.
- Check the device and all its accessories before each use.
- The device should only be used under the constant supervision of a trained and authorized practitioner.
- You must have been trained by an authorized person from SWIMS S.A.S./WINBACK before using the device.
- Keep the device out of the reach of children.
- Improper use of the controls or settings, or the carrying out of procedures other than those specified herein, may result in damage to the equipment or its accessories.
- Choose the position of the neutral electrode carefully. Never place the return plate on an area of the patient's skin that is delicate or in a bony area.
- The values displayed on the main screen may vary between +/-10%.
- The emergency stop button must always be plugged in and accessible to the patient.
- Always make sure that the patient does not have any contraindication before using the device.
- To use hands-free mode with fixed electrodes, take care not to exceed an intensity of 40% and use LOW mode.

Technical:


- The device and its accessories must be checked regularly to ensure that they are functioning properly. It is recommended to have the equipment serviced by the manufacturer every two years.
- The device must be used exclusively with the accessories provided. Using inappropriate accessories might lead to electrical burns or cause the device to malfunction.
- The device must be used under conditions of standard atmospheric pressure (around 1013 hPa), within a temperature range of -10°C and 40°C, and within a relative humidity range of 20% to 90%.
- To avoid the risk of electric shock, always unplug the device from the electrical outlet before carrying out any maintenance procedures*.
- Particular care should be taken when using the device in proximity to other equipment that may generate electromagnetic interference. Avoid using it around other electrical equipment.
- Never use accessories that are cracked or broken.
- Ensure that your electrical network complies with the regulations in force.
- Never bend the power cord or the cables of the accessories excessively.
- Keep liquids away from the device.
- Do not set up the device in a place where it is difficult to remove the power plug.
- A stable wifi connection is mandatory

***WARNING:** Only SWIMS S.A.S./WINBACK technicians and/or technicians approved by the company are permitted to open the device. Users are strictly forbidden to open the device or try to change any components..

3. DEVICE DESCRIPTION

3.1 LOCATION OF BACK4 CONTROLS

Note: Dimensions are in millimeters

3.2 MAIN SCREEN/FUNCTIONS

In order to provide treatment that is most appropriate to the patient, BACK4 offers tecartherapy, electrostimulation, or a combination of the two.

- FR
- EN
- ES
- IT





Item	Function
Mode selection	Select the mode, the frequency and the signal modulation*.
Signal intensity	Set the intensity level of the output signal.
Settings	Access to parameters (serial number, software version, WIFI configuration, sound configuration).
Wifi	WIFI connection. When the device is connected to WIFI, the signal turns green. After 5 minutes of inactivity, the WIFI turns off and grey.
Help	This explains the different effects related to each mode.
Treatment intervention buttons	Start, pause and stop the treatment.
Time selection	Select the total treatment time. Displays the total treatment time and the time remaining.
Time mode	Total time increment for each mode.

*Note: For details of each mode, please refer to section 4.2, Modes and functions

3.3 BACK4 ENERGY

The device can be used differently depending on the needs of the patient and the recommendations of the professional.

Mode	Frequency	Diagram																						
RET	RET: monopolar resistive energy at 500kHz RET + Hi-TENS : monopolar resistive energy at 300kHz																							
CET	Soft: monopolar capacitive energy at 500kHz Deep: monopolar capacitive energy at 300kHz																							
MIX	MIX: multipolar resistive energy via a sinusoidal wave at 1000kHz MIX + Hi-TENS: multipolar resistive energy via a sinusoidal wave at 500kHz	 <table border="1" data-bbox="721 1079 1447 1182"> <thead> <tr> <th>Cycle</th> <th>10%</th> <th>20%</th> <th>30%</th> <th>40%</th> <th>50%</th> <th>60%</th> <th>70%</th> <th>80%</th> <th>90%</th> <th>100%</th> </tr> </thead> <tbody> <tr> <td>ms</td> <td>20</td> <td>40</td> <td>60</td> <td>80</td> <td>100</td> <td>120</td> <td>140</td> <td>160</td> <td>180</td> <td>200</td> </tr> </tbody> </table>	Cycle	10%	20%	30%	40%	50%	60%	70%	80%	90%	100%	ms	20	40	60	80	100	120	140	160	180	200
Cycle	10%	20%	30%	40%	50%	60%	70%	80%	90%	100%														
ms	20	40	60	80	100	120	140	160	180	200														
Hi-TENS	The pulse frequency can be set at 2, 5 or 25 Hz. The carrier frequency of each pulse is 300 kHz for Hi-TENS (RET) mode and 500 kHz for Hi-TENS (MIX) mode.																							
Hi-EMS	Electrical stimulation with modulated frequencies of 1.5kHz and 4kHz.																							

All diagrams were obtained using a 500 ohm load, except for MIX mode which uses a 50 ohm load.

3.4 ACCESSORIES

Your device will be delivered with a range of accessories. All these accessories are meant for external use. Never use your device with an accessory which is not on the following list.

FR

EN

ES

IT

Tecar X handset (x2)
B4_ACC.X

RET Bracelet (x2)
(applied part)
REF: B4_ACC.34

Neutral return plate
(applied part)
REF: B4_ACC.1

Neutral handle
(applied part)
REF: B4_ACC.41

RET electrodes (applied part)

Ø 70
REF: B4_ACC.29

Ø 60
REF: B4_ACC.28

Ø 40
REF: B4_ACC.27

Convexe 40
REF: B4_ACC.22

CET electrodes (applied part)

Ø 40
REF: B4_ACC.23

Ø 60
REF: B4_ACC.24

Ø 70
REF: B4_ACC.25

Convexe 60
REF: B4_ACC.30

Hi-RET electrode
(applied part)
REF: B4_ACC.45

Multipolar electrode S
(applied part)
REF: B4_ACC.36

Multipolar electrode L
(applied part)
REF: B4_ACC.37

RET Fixed pad
(applied part)
B4_ACC.13

Neutral Fixed pad
(applied part)
REF : B4_ACC.14

TECAR mobile handle
REF : B4_ACC.12

RET adhesive cable
B4_ACC.15

Neutral adhesive cable
B4_ACC.16

RET/CET generic cable (x2)
REF : B4_ACC.7

Neutral generic cable
REF : B4_ACC.9

Emergency stop button
REF : B4_ACC.2

Power cord
REF : B4_TEK.1

Winback trolley
REF: TROLL_NEW

Note: Generic neutral and RET/CET cables are for use with the corresponding accessories. It is mandatory to connect one side of the cable to the accessory and the other side to the device.

CONSUMABLES

**Antiseptic wipes
(1 pack of 30 units)**
REF : CONS3.0

Winback towel
REF : MKT8.0

**Essential conductive
cream (1L)**
REF : CONS1.2

**Return adhesive plate
(1 pack of 5 units)**
(applied part)
REF : CONS2.0

4. OPERATING PRINCIPLES

4.1 EFFECTS OF THE HIGH-FREQUENCY CURRENT

BACK4 generates high-frequency sinusoidal current with monopolar and multipolar modes. It can be used with one or two active accessories and one neutral accessory. When the two/three accessories are in contact with the patient, the electrical circuit closes and current is delivered.

The effects of radiofrequency (RF) currents are related to their frequency and intensity.

All living tissues have some degree of resistance or impedance to the flow of current. Scientific studies have shown that the frequency band between 300-1000kHz induces an increase in intra- and extracellular exchanges.

BACK4 has two main modes: CET (Capacitive Electrical Transfer) and RET (Resistive Electrical Transfer). CET mode uses less energy than the RET mode (the conversion of electrical energy to increased internal temperature is much higher in RET mode). CET mode has a more superficial action and targets soft tissues, while RET mode has a deeper action and targets hard tissues. BACK4 also has a MIX mode (multipolar) for superficial action on the epidermis and dermis. BACK4 can be used with an electrical stimulation for muscle stimulation and drainage.

4.2 MODES AND FUNCTIONS

Capacitive mode (CET)

Monopolar capacitive energy at 500kHz (Soft) or 300kHz (Deep). This mode is used for superficial action on soft tissues with a high water density (muscles, fascia, viscera, lymphatic, vascular, etc.). You can also select Dynamic mode which alternates the frequencies every 4 seconds.

Resistive mode (RET)

Monopolar resistive energy at 500kHz. The responses produced in resistive mode are concentrated in the deepest and most resistant tissues, such as the joints and tendons. These areas of high density and low water content focus more energy than in areas with greater water content. RET mode creates an internal thermal increase based on the direct application of the current on the skin (because the electrode is not isolated). This action allows a greater temperature increase to be achieved using less energy. To use RET at 300kHz, select Hi-TENS mode with RET mode.

Hi-TENS mode (RET)

This mode allows the use of tecartherapy with pulsed current at the maximum intensity level. The Hi-TENS frequency can be set to 2, 5 or 25Hz.

This mode will be used for an analgesic effect, relieving pain and releasing tension. The 2Hz and 5Hz frequencies can be used to relieve chronic pain, while the 25Hz frequency can be used to treat acute pain. “Dynamic” mode, which alternates frequencies every 4 seconds, can also be selected.

Hi-EMS mode

The frequencies used in this mode are 1.5 kHz focal and 4 kHz radial (modulated frequencies).

This mode has an excito-motor action.

This medium frequency current allows a better recruitment of the motor units with a deep action.

Multipolar mode (MIX)

Multipolar resistive energy at 1,000kHz. The reactions produced by MIX mode mainly focus on the tissues that are in contact with the electrode because this mode operates across a very short circuit, with both the mobile and fixed electrodes contained in the MIX electrode. The action is superficial, affecting the epidermis and dermis, and suitable for cutaneous and subcutaneous treatments. To use MIX at 500kHz, select Hi-TENS mode with MIX mode.

Hi-TENS mode (MIX)

This mode can only be used in combination with MIX mode. It allows the use of tecartherapy with pulsed current at the maximum intensity level. The pulse frequency can be set to 2, 5 or 25Hz.

This mode will be used for an analgesic effect, relieving pain and releasing tension. The 2Hz and 5Hz frequencies can be used to relieve chronic pain, while the 25Hz frequency can be used to treat acute pain. “Dynamic” mode, which alternates frequencies every 4 seconds, can also be selected.

Signal modulation

“Low” mode: The signal is active for 0.5 seconds, then inactive for 0.5 seconds to enable the practitioner to work athermally. Low mode permits the maintenance of a specific biological stimulus while limiting diathermy, which is particularly useful during acute inflammation phases. This is made possible by the controlled alternation between emission of energy targeting the tissues followed by a neutral period which promotes the physiological maintenance of the local temperature.

“Medium” mode: TECAR signal at the selected intensity without modulation.

“Boost” mode: The signal is active for 0.5 seconds at the selected intensity, followed by a further 0.5 seconds at +30% of the selected intensity. The principle of this mode is to increase the intensity for a very short and repetitive period. The effect is to force the poorly vascularized areas to get back into motion. Feedback from users attests to a significant time-saving. It alleviates tension and releases fibroses twice as fast.

To select these modes, press the ‘Mode’ button. Options will appear on the right side of the interface. Select the desired mode.

“SWAP” mode : Alternation between RET mode and Hi-TENS mode for 2 seconds. Can only be used with RET and Hi-TENS modes.

Additional options

Selection of intensity and treatment time: If you press on the intensity or the treatment time, a numeric keypad appears on the right of the screen. Select the intensity (from 1 to 100%) or the time (from 1 to 60 minutes) then press to the validation button.

In double RET mode or in RET+CET:

- If you are using a neutral accessory in return, select the “Neutral Fix” mode.
- If you are using a return RET electrode, select the “Neutral Mobile” mode.

5. PERFORMING A TREATMENT

5.1 PREPARING THE DEVICE AND PERFORMING A TREATMENT

After making sure that the patient does not have any contraindications and they fully understood the purpose of the treatment, choose the appropriate protocol for the pathology being treated and make sure that the device and its accessories are in perfect working order.

Step one: Clean the accessories and the area(s) to be treated. The area should be clean and dry.

Step two: Position the device so that the air vents at the back and sides are not obstructed. The cooling system is air based and it is vital that the vents are not obstructed. Place the device on a flat surface.

It is now ready to be used.

5.2 TREATMENT

1. Rear panel and connections

For details of the different connectors, please refer to the rear view of the device, chapter « 3.1 Location of BACK4 controls ».

Turn on the device (rear power button + front ON/OFF button).

Connect the accessories (from 1 to 3 as needed) and the emergency button and check if the connectors of the accessories are correctly connected to the device.

FR

EN

ES

IT

2. Set treatment parameters

Before carrying out a treatment, set the parameters based on the patient and the pathology. Select the mode, the signal intensity level and the duration. When an accessory with an electrode is connected, the default mode corresponding to that electrode is automatically activated, and only the modes that it supports can be selected:

- When a mode is available, an illuminated line appears under the mode.
- When a mode is not available, this line is not illuminated and the entire column is grayed out.
- When a mode is available but not activated, it is not highlighted.
- When a mode is activated, the entire column is highlighted.

In the image below, only CET, RET, Hi-TENS and Hi-EMS modes can be selected (illuminated line), while CET, RET and Hi-TENS modes are activated (columns are highlighted). MIX and Hi-TENS modes cannot be activated (modes are grayed out).

With the TECAR X accessory, you can navigate on the interface. The TECAR X accessory logo appears on the selectable mode:

- To change mode, swipe twice to the left or to the right.

- To activate/deactivate the mode, make a pronosupination movement with the wrist.

- To change the intensity of the treatment, tap the accessory with your finger. To increase the intensity by 10%, tap 3 times, and to decrease the intensity by 10%, tap twice.

To use this function, please activate the "TIC" in the "Settings" menu.

- In double RET or CET+RET, “shock” the two electrodes together to switch to “Neutral Fix” or “Neutral Mobile” mode.

FR

EN

ES

IT

3. Conductive cream


A conductive cream or gel (for mucous membranes) must be used to facilitate electrical exchanges with the accessory and ensure a sufficient level of safety for the patient's skin. WINBACK has developed a cream and gel to meet these requirements. If the cream/gel used during treatment does not meet certain minimum conditions of electrical conductivity, the results obtained will not be optimal and the risk of side effects will be increased.

Warning: This technology requires a specially formulated cream/gel. A cream or gel that has not been formulated for use with the equipment will form a crown around the electrode as the temperature increases. Conversely, a suitable product will spread evenly over the skin under the same conditions. Only Winback conductive cream should be used.

4. Connection/disconnection of electrodes

The TECAR X and TECAR mobile accessories can be adapted with the addition of an electrode. To connect the two, clip the electrode on to the handle where the mark indicates. To disconnect them, simply pull on the electrode.

Caution: Electrodes must be clipped in place and are not meant to be screwed or unscrewed.

5. Position of the electrodes

In some cases, a loss of contact between the two electrodes can occur because the distance is too great. The electrodes should be placed according to the area to be treated.

The position of the patient depends on the area being treated. It is used to determine the best position for the return accessory and the active accessory.

Example: for treating the knee, the neutral accessory should be placed under the thigh

6. Application of the electrodes

Sufficient pressure must be exerted to maintain optimal electrical contact over the entire surface of the electrode, while allowing smooth movement.

The electrode should be moved in small circles using a slow, uniform motion, and without stopping. It is preferable to reduce the intensity if it is uncomfortable for the patient, then increase it again gradually, rather than maintaining a high power and being forced to move the electrode faster.

When the electrode is not in contact with the patient's skin, the timer stops, and no current passes through the electrode.

Warning: The surface of the electrode must be in perfect contact with the skin. The electrode should not be tilted, but instead held completely flat during treatment to avoid any discomfort related to overheating, the electrical arc or burns. Always make sure that your patient does not feel any sensation or pain that is unusual.

Insufficient contact

Good contact

7. Additional functions during treatment

When using the TECAR X, the color of the mode selected matches the color of the illuminated logo on the accessory:

- The logo flashes when treatment is paused
- The logo pulses during treatment
- The logo remains illuminated but does not change in brightness when the electrode is not in contact with the skin or when the treatment is stopped

8. Connected accessories and possible modes of use

	CONNECTED ACCESSORIES		
	CET accessory	RET accessory	Multipolar accessory
OPTIONS FOR USE	CET	RET Hi-TENS (RET) Hi-EMS RET + Hi-TENS (RET) RET + Hi-EMS Hi-TENS (RET) + Hi-EMS RET + Hi-TENS (RET) + Hi-EMS	MIX MIX + Hi-TENS (MIX)

6. MAINTENANCE OF THE DEVICE

6.1 GENERAL MAINTENANCE

To avoid any risk of electric shock, the device must be turned off and unplugged before cleaning or technical maintenance is undertaken.

Never use corrosive products or an abrasive cloth during cleaning, as these may damage the device. Damage caused by the use of products other than those that are recommended, or resulting from inappropriate cleaning, is not covered by the warranty.

Never pour liquids directly on the screen or the device. Do not clean the screen with detergents (the presence of liquid could cause a malfunction). If cream or liquid is spilled on the device, dry it immediately to avoid damage.

The BACK4 and its accessories are reusable. To avoid potential contamination of patients by pathogens, it is essential that the device and its accessories are clean.

The following cleaning protocol should be carried out after each treatment:

- Wipe any excess cream from the accessories.
- Carefully clean the accessories using an appropriate disposable disinfectant wipe*.
- Dry the accessories using a clean, disposable paper towel.
- Make sure that all accessories are clean, free from conductive cream and completely dry.

*Always ensure that the disinfectant you use is suitable both for a medical device and for skin contact (e.g.: appropriate products from Bactynea, Stericid, Purell or Aseptonet).

Part	Cleaning material(s)	Cleaning process
Touchscreen	A clean, soft cloth lightly moistened with rubbing alcohol (surgical spirit).	Clean the touchscreen with the cloth, allowing any excess alcohol to dry naturally.
Power cord	A clean, dry cloth.	Wipe the power cord regularly with a dry cloth. Moisture and dust can cause fire or electric shock.
Accessories and electrodes	A clean, soft cloth lightly moistened with rubbing alcohol (surgical spirit).	After each treatment, wipe the accessories with a soft, dry cloth or a slightly damp microfiber cloth, then with a disinfectant wipe.

Note: If you have not used the device for a few days, apply the same cleaning protocol before using the BACK4 again.


WARNING: Repair or technical maintenance of the device, or any of its accessories, are prohibited. Only SWIMS S.A.S./WINBACK technicians are permitted to open or service the device. Users are strictly forbidden to open the device.

6.2 REMOTE MAINTENANCE

This device is equipped with a remote maintenance module to allow easy access to software updates via wifi. The connection is only used for device updates and no personal information is collected.

When the device is switched on for the first time, it will offer to connect to wifi. Select a wifi network, enter the password and validate.

Once connected to wifi, press X to return to the settings screen.

You can check the WIFI connection and signal strength using the WIFI logo. When the device is connected to wifi, the symbol turns green.

When a software update is available it will be visible in the settings tab, with the option to proceed or not. To proceed with the update, follow the instructions on the screen.

Note: the WIFI automatically disconnects after 5 minutes of inactivity. No personal data is collected.

7. TROUBLESHOOTING & ERROR MESSAGES

7.1 TROUBLESHOOTING

WARNING: Never attempt to open the device or use tools on the device or its accessories. If the appropriate troubleshooting procedure does not resolve the issue, contact the service team.

- Everything is set up ready for treatment. On pressing START, the button is active but the timer does not start counting down.

Check whether the issue has been resolved after completing each of the following steps.

Plug in and unplug all connectors and make sure that none of them are damaged.	If any of the cables/connectors are damaged, call after-sales service for replacements.
Swap the neutral accessories.	If this works, call after-sales to request a replacement neutral electrode.
Add a little more cream to the treatment zone.	The conductive cream is essential for insuring sufficient contact between the accessories and the patient's skin.

7.2 ERROR MESSAGES

Error	Problem	Cause	Solution
254	Emergency button pressed	<ol style="list-style-type: none"> 1. The patient presses on the button. 2. The button is misconnected. 3. The button is damaged. 	<ol style="list-style-type: none"> 1. Make sure that the patient did not press the button OR check the patient's skin and make sure that everything is ok. Turn off and restart the device using the button on the rear. 2. Unplug the button and plug it in correctly. If the problem persists, call after-sales service. 3. Make sure that the button is not damaged. If it is, call after-sales service.
115	Neutral Fix mode selected	The "Neutral Fix" mode can only be used with a neutral accessory in return.	Use a neutral accessory in return when "Neutral Fix" mode is selected, or select "Neutral Mobile" mode.
119	MIX accessory connected	Another accessory is connected in addition to the MIX accessory.	Either disconnect the other accessory to use the MIX accessory, or disconnect the MIX accessory to use the other accessory.
218	Boards version not compatible on way 1 or 2	The board software version is not similar to the GMU version.	<ol style="list-style-type: none"> 1. Please update the software version by pressing the "UPDATE" button in the "Settings" menu. 2. Check that the version of TECARX accessories is the same as the version of BACK4. 3. If problem persists, call after-sales service.

For all other problems or errors, please turn off the device and contact after-sales service.

WARNING: You should never open the device.

8. TECHNICAL DATA

8.1 TECHNICAL SPECIFICATIONS

General	
Power input	100-240V a.c. 50/60Hz, 300VA
Max. power consumption	300W
Fuse	2 x (2A, 250Ac 50x20)
Classification (EN 60601-1)	II
Size	403.6x379.0x190.5mm
Weight	6kg
Medical device class	Ila

Output parameters	
Frequency	1,5kHz, 4kHz, 300kHz, 500kHz, 1MHz
Treatment timer range	1 to 60 minutes
Safety	Zero start, safety switch

FR

EN

ES


IT

Environmental conditions for use	
Temperature	-10°C to 40°C
Humidity	20% to 90%
Pressure	800 to 1060hPa

Environmental conditions for storage and transport	
Temperature	-20°C to 70°C
Humidity	10% to 95%
Pressure	800 to 1060hPa

Wifi	
Frequency band	2400-2483.5MHz
Max. power transmitted	Up to 20dBm

Label and packaging

Information about symbols

- 
 Manufacturer: WINBACK GO EAST - #608, 200, Geopdosi-ro, Jijeong-myeon, Wonju-si Gangwon-do, South Korea.
- 
 Date of manufacture
- 
 Serial number
- 
 It is mandatory to read the user manual before using the device
- 
 Type BF applied parts
- 
 The device must be returned to SWIMS S.A.S./WINBACK at the end of its life in accordance with the DEEE 2002/96/CE and EN50419 European standards
- 
 SWIMS SAS is the legal representative of WINBACK GO EAST in the European Community
- 
 Do not sit on the device
- 
 Do not step or climb on the device
- 
 Device generates non-ionizing emissions
- 
 This way up. Do not bend the packaging
- 
 Contents are fragile and must be handled with care
- 
 Keep away from water
- 
 Handle with care
- 
 Temperature limits
- 
 Humidity limits
- 
 Atmospheric pressure limits
- 
 Warning: electrical output

8.2 COMPOSITION CONFORMITY STATEMENT

The BACK4 device does not incorporate, as an integral part, any medicinal product, human blood or derivatives. The device does not incorporate tissue of animal origin.

8.3 TECHNICAL MAINTENANCE

If any parts or accessories are defective, please contact SWIMS S.A.S./WINBACK after-sales service (sav@swims.team).

The life expectancy of the BACK4 device is 5 years. However, to guarantee maximum efficiency, maintenance must be carried out every 2 years.

During this maintenance, output power and general performance will be checked to ensure the efficiency and safety of the device are still compliant with the manufacturer's standards.

9. END OF LIFE

Reflecting prevailing regulations, SWIMS S.A.S./WINBACK supports the DEEE Pro ecosystem recycling program, which collects used medical material free of charge.

At the end of its life, return your BACK4 to: SWIMS S.A.S./WINBACK - Les Bureaux du Parc - 2474 RD 6007 - 06270 Villeneuve-Loubet, France. If you are located outside Europe, contact SWIMS S.A.S company - for details of the address to return your device to.

10. SAFETY AND ELECTRICAL STANDARDS

The BACK4 device is a Class IIa medical device with Class C software.

To ensure optimal safety and efficiency, the BACK4 has been developed and manufactured in accordance with the following standards:

Name	Description
IEC 60601-1	General requirements for the safety of medical electrical systems
IEC 60601-1-2	Electromagnetic compatibility for medical devices
IEC 60601-2-10	Medical electrical equipment - Part 2-10: specific requirements for basic safety and essential performance of nerve and muscle stimulators
ISO 14971:2014	Application of risk management to medical devices
93/42/EEC	Essential requirements for medical devices
ISO 13485:2016	Quality management systems for medical devices
2014/53/EU	Radio equipment directive (RED)

In accordance with Directive 93/42/EEC, the following tables provide information regarding the EMC characteristics of this electromedical equipment.

For details of all the tests carried out on BACK4, refer to the electrical safety and EMC testing report with regard to ISO 60601-1 and ISO 60601-1-2.

Guidance and manufacturer's declaration - electromagnetic emissions		
The BACK4 is intended for use in an electromagnetic environment with the characteristics set out below. The user must ensure that it is used in such an environment.		
Emissions	Compliance	Electromagnetic environment guidance
RF emissions CISPR11	Group 2	Its RF emissions are very low and are not likely to cause interference with nearby electronic equipment.
RF emissions CISPR11	Class A	The electromagnetic equipment is suitable for use in all settings, other than domestic settings and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000 -3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity			
The BACK4 is intended for use in an electromagnetic environment with the characteristics set out below. The user must ensure that it is used in such an environment.			
Immunity test	Test condition	IEC 60601 Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, cement, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
RATED power frequency magnetic fields IEC 61000-4-8	50Hz or 60Hz	30 A/m (rms)	Power frequency magnetic fields should be equivalent to those in a commercial or hospital setting.
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV @ 10kHz Repetition frequency	± 2 kV @10kHz for power supply lines	Mains power quality should be that of a typical commercial and/or hospital environment.
Surge IEC 61000-4-5	+/- 0.5kV, 1kV/L-N +/- 0.5, 1, 2kV/ L-PE +/- 0.5, 1, 2kV/ N-PE	+/- 0.5kV, 1kV/L-N +/- 0.5, 1, 2kV/ L-PE +/- 0.5, 1, 2kV/ N-PE	Mains power quality should be that of a typical commercial and/or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% UT; 0° 0% UT; 70% 0% UT; 0%	0,5 cycles 1 cycle 25/30 Cycles (50/60Hz) 250/300 Cycles (50/60Hz) (5s)	Mains power quality should be that of a typical commercial and/or hospital environment. If the user of the BACK4 requires continued operation during mains power interruptions, it is recommended that the device is powered from an uninterruptible power supply.

Conducted RF induced by RF fields IEC 61000-4-6	3Vrms outside the ISM band 6Vrms in the ISM band	3Vrms outside the ISM band 6Vrms in the ISM band	Mains power quality should be that of a professional healthcare facility environment.
<p>Comment:</p> <p>a) Not applicable to CLASS II medical electrical equipment and medical electrical systems.</p> <p>b) The ISM (industrial, scientific and medical) bands between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40,66 MHz to 40.70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7,3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.</p>			

Guidance and manufacturer's declaration - electromagnetic immunity

The BACK4 is suitable for use in the specified electromagnetic environment. The customer and/or the user of the BACK4 should assure that it is used in an electromagnetic environment as described below.

Immunity test	IEC 60601-1-2 Test Level	Compliance Level	Electromagnetic environment guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3Vrms outside the ISM band 6Vrms in the ISM band</p> <p>3V/m (80MHz-2.7GHz)</p>	<p>3Vrms outside the ISM band 6Vrms in the ISM band</p> <p>3V/m (80MHz-2.7GHz)</p>	<p>Portable and mobile RF equipment should not be used at a distance less than the recommended separation distance, calculated using the equation for the transmitting frequency.</p> <p>Recommended separation distance:</p> $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2.3 \sqrt{P} \text{ 800 MHz to 2,7 GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey(a), should be less than the compliance level in each frequency range(b).</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

- a. Field strengths from fixed transmitters, such as base stations for radio telephones (cellular/cordless) and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be accurately predicted theoretically. To assess the electromagnetic environment resulting from fixed transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment and/or system is being used exceeds the applicable compliance level above, the BACK4 should be checked to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the BACK4.
- b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

NOTE : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Recommended separation distances between portable and mobile RF communications equipment and the above listed models

The BACK4 is intended for use in an electromagnetic environment in which radiated RF disturbance variables are controlled. The customer or user of the BACK4 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BACK4 as recommended below, based on the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Protection distance based on transmitter frequency (m)		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: From 80 MHz to 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

WARNING:

Use of this equipment next to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be monitored to verify that it is operating normally.

Portable RF communications equipment (including peripherals such as antenna cables and external antennae) should be used no closer than 30cm (12 inches) to any part of the electromagnetic equipment or system, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

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SWIMS SAS - WINBACK